

DSC using Q20 instrument

Experimental Considerations

- 1) Before starting the instrument consider your experiment conditions:
 - a. What temperature range are you testing?
 - i. Maximum temperature on this system is 550 °C
 - ii. Minimum temperature is -90 °C
 - b. What scan rate will you be using, will you need different scan rates in different temperature zones?
 - i. Depending on your goal temperature there is a maximum cooling rate that can be achieved
 - ii. From 550 °C to 300 °C; 100 °C/min
 - iii. From 550 °C to 120 °C; 50 °C/min
 - iv. From 550 °C to -20 °C; 20 °C/min
 - v. From 550 °C to -50 °C; 10 °C/min
 - vi. From 550 °C to -75 °C; 5 °C/min
 - vii. From 550 °C to -90 °C; 2 °C/min
 - c. Will your sample react with the sample pan at the desired temperatures?
 - i. If sample reacts with pan this will damage the instrument, **stop here**
 - ii. If you're not sure then **stop** and go look up more information
 - d. How will your sample break down under heating? Will it decompose to form gaseous by products?
 - i. Samples that off gas when heated require special pans and setup
 - ii. If you are not sure how the sample will react under heating test on the TGA to determine decomposition temperatures and when weight losses occur is recommended
 - iii. If in doubt go to the literature to see how others have handled samples of similar type

Sample loading

- 1) Select a pan and a lid for your sample
- 2) Make sure that you have a reference pan of the **exact** same pan type
- 3) Place the pan and lid on the blue colored pan holder




- 4) Put pan/lid/holder setup on analytical balance and zero the balance



- 5) Carefully load the pan with sample
 - a. The sample should lay flat on the bottom of the pan
 - b. Make sure there is no sample outside of the inside part of the pan before weighing
- 6) Check weight of the sample
 - a. Mass should be between 1 mg and 20 mg
- 7) Once appropriate amount of sample is loaded carefully place lid on pan so that it sits level



- 8) Place the pan/lid/pan holder setup in the pan crimper and pull the lever down completely then raise it
 - a. The whole setup should sit snugly in place before you pull the lever




- 9) Check the pan and lid to see if the lid is fully seated in the pan with no visible sample
- 10) Weigh the sample after sealing the pan, this is the mass that will be input in the instrument program
- 11) Load sample pan and reference pan into the sample cell in the instrument
 - a. To load pans remove heat cover, shown below



- b. Next remove both cell lids using tweezers to handle them




- c. The sample pan should be placed on the column closest to you while the reference pan should be be on the column towards the back



- d. Replace all of the cell lids and heat cover before setting up computer system and starting the scan

Computer Program Setup

- 1) Open the N₂ tank attached to the DSC
- 2) On the instrument computer open the instrument control program by clicking on the instrument icon outlined in red in the image below



- 3) In the instrument program go to **Control** select **Event** and click **On** to turn on the cooling system
 - a. Do not start any experiment until the Flange temperature reaches below $-70\text{ }^{\circ}\text{C}$
 - b. Flange temperature is displayed in the table on the right hand side of the screen


Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	41.31 °C
Heat Flow	-0.016 mW
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	41.53 °C
Heater Temperature	41.64 °C
Sample purge flow	0.00 mL/min

Running Segment Description

- 1 Sampling interval: 0.10 s/pt
- 2 Equilibrate at -50.00 °C
- 3 Ramp 10.000 °C/min to 40.00 °C

Heat Flow (mW) vs Temperature (°C) graph showing a flat line at approximately 0.00 mW.

- 4) While the flange is decreasing set up experimental conditions
- 5) Go to the Summary Tab
 - a. Under **Mode** select **Standard**
 - b. Under **Test** select **Custom**
 - c. Under **Sample Information** fill in the sample name, select the pan type, input sample mass in mg from before and add any comments in the comments section
 - i. If you are using the supplied DSC pans select Aluminum
 - d. Next to **Data File Name** click on the button with the box on it to open the file save menu
 - e. In this menu select the folder you want to save in and fill in the file name then click save



- 6) Go to the Notes tab
 - a. In the operator field type in your name
 - b. In the purge gas section select the gas you will be using (N₂ or air)
 - c. The Flow rate should be kept at 50 mL/min

QSeries - [Q20-2870 - DSC Q20@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Standby Temp: 41.31 °C Sig 1: 41.31 °C Sig 2: -0.016 mW Sig 3: 0 min Store on Event on

Experiment
Standard Sequence

Sequence No. 73
Run 1:

Summary Procedure Notes

Notes

Operator: Jen

Extended Text

Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Auto Analysis

Autoanalyze

Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	41.31 °C
Heat Flow	-0.016 mW
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	41.53 °C
Heater Temperature	41.64 °C
Sample Purge Flow	0.00 mL/min

#	Running Segment Description
1	Sampling interval: 0.10 s/pt
2	Equilibrate at -50.00 °C
3	Ramp 10.000 °C/min to 40.00 °C

Heat Flow (mW)


Temperature (°C)

01 15:40 min Append Apply Cancel Help


Ready Stand by Standard Seg 0 in Run 1 12:13:27

7) Go to the Procedure tab

- a. In the drop down next to **Test** select **Custom**
- b. In the notes section you can fill in any information about your test that you would like to have
- c. Under **Method** click on the **Editor** button to change the method




d. This will open the method editor window shown below



- e. Click on the **New Method** button (White sheet icon) to start a new method
- f. Click on and drag the segment types you want to use in your experiment into the white box on the left
- g. Most common segments used are:

- i. **Ramp**: Changes furnace temperature at a given rate in °C/min to a specified temperature
- ii. **Isothermal**: Holds furnace at last set temperature for a specified amount of time in minutes
- iii. **Equilibrate**: Heats or cools the system as quickly as possible to a particular temperature and then holds until the temperature stabilizes

8) Once you have confirmed that all method and data file information has been properly input click on **Apply** in the bottom left side of the screen



Starting scan

- 1) Check that sample pan and reference pan are properly loaded in the instrument cell
- 2) Check that the desired method is properly loaded and does not exceed the max temperature of 550 °C
- 3) Check that the Flange Temperature is below -70 °C
- 4) If all of these are set up then start the run by clicking the green play button

Post Experiment tear down

- 1) Once then scan is complete the sample pan can be removed
 - a. If you are doing multiple samples then replace the finished sample with the next sample and start that scan
 - b. If you are done with your experiments then just remove your sample pan and replace the cell covers and heat shield
- 2) When you are done with the instrument turn off the cooling system
 - a. Got to **Control** then **Event** and click **Off**
 - b. You will hear the cooling system shut off
- 3) Leave the N₂ gas tank open until the Flange Temperature warms up to higher than 10 °C
 - a. This is to prevent forming ice in the sample cell
 - b. If ice forms in the cell then it will have to be baked out and recalibrated resulting in the instrument being down for several days