

# E-Cylinders in the Hospital

## Full Tanks

- Will continue to be delivered to the tank room
- DOWL 0140



## Partially Full

- As soon as tanks or cylinders leave the tank room they are considered "**Partially Full**"



## In Use

- Tanks deemed "**in use**" when on a stretcher, bed, wheelchair, or code cart



## Empty

- When a cylinder is "empty" (below 500 psi), they must be placed in a carrier.
- Carriers are located on each floor in the room by the elevators (ex: 4206, 5206, 6206)
- 3rd Floor "empty" carrier is located in 3145.

## Filled Carrier

- When carrier is filled with "**Empty**" cylinders, it **MUST** be taken to the tank room on Dowling Hall Dock, Room 0140
- Place "**Empty**" cylinders in appropriate location in tank room

## Labeling

- IMPORTANT:** Do not co-mingle "Empty Cylinders" and "Partially Full Cylinders" on the floors
- Each cylinder holder is properly labeled as such
- Transportation carts should be labeled with a "Partially Full" sticker. Do not leave empty cylinders in carts.

