


Name of Policy:	<u>LTME/Video EEG monitoring checks</u>	
Policy Number:	3364-138-20	
Department:	Neurodiagnostic Services	
Approving Officer:	Senior Hospital Administrator	
Responsible Agent:	Director, Pulmonary Services	
Scope:	The University of Toledo Medical Center Neurodiagnostic Services	
		Effective Date: 6/1/2023 Initial Effective Date: 5/27/2014
<input type="checkbox"/> New policy proposal <input type="checkbox"/> Major revision of existing policy		<input checked="" type="checkbox"/> Minor/technical revision of existing policy <input type="checkbox"/> Reaffirmation of existing policy

Policy Statement

To establish guidelines for checking running LTME/Video EEG monitoring patients.

Procedure

All patients on LTME/Video EEG monitoring will be checked at least twice daily for the following.

1. Impedance
2. EKG patches.
3. Check the connection of imaging friendly electrodes.
4. Check for active connection to the server, if available.
5. Make sure EEG is recording.
6. The camera should be position to appropriately visualize patient.
7. Check for any sedation running.
8. Check that cables and cords are not tangled.

Document on the recording:

1. Impedance check.
2. That cable is correctly connected to the breakout or EMU40 box.
3. Sedatives used and dosage.

Approved by:		Review/Revision Date:
<u>/s/</u> Michael Taylor Director, Pulmonary Services	<u>6/22/2023</u> Date	05/27/2014 05/29/2020
<u>/s/</u> Russell Smith Senior Hospital Administrator	<u>6/29/2023</u> Date	
		Next Review Date: 6/1/2026
Policies Superseded by This Policy:		